

Step 1 Fill In Your Information

Dentist Name: _____ Office Name: _____
 Patient Name: _____ Address: _____
 City: _____ State: _____ Zip: _____ Phone: _____
 E-Mail: _____ Fax: _____

Step 2 Need a Phone Call?

Yes No (phone call may delay production time)

Step 3 Select a Device

Refer to website catalog for additional information on a specific device.

- Respire Blue+** (Dorsal Fin)
 - Hard Hard w/ Reduced Lingual Coverage Hard/Soft
- Respire Blue EF+** (Endurance Framework)
 - Hard w/ Ball Clasp Hard w/ Hard with C Clasp
- Respire Pink** (Herbst)
 - Hard Hard w/ Reduced Lingual Coverage Hard/Soft
- Respire Pink Micro**
 - Hard Hard w/ Reduced Lingual Coverage Hard/Soft
- Respire Pink EF** (Endurance Framework)
 - Hard w/ Ball Clasp Hard w/ C Clasp
- Respire Pink EF Micro** (Endurance Framework)
 - Hard w/ Ball Clasp Hard w/ C Clasp
- TMD Treatment**
 - Gelb Lingual Coverage
 - Full Coverage Upper Full Coverage Lower
 - NightGaurd Upper Hard NightGaurd Upper Hard/Soft
 - NightGaurd Lower Hard NightGaurd Lower Hard/Soft
 - Other (Step 7: Special Instructions)

Step 4 Select Add-ons (Only for Sleep Appliances)

- Anterior Bump (Discluding Element) \$25 Elastic Hooks \$25
- DentiTrac® Compliance Chip (Respire Pink (Herbst) Devices Only) \$120
- Custom Bite Repositioner \$89

Step 5 Check the Enclosed Items

Impression: Upper Lower Bite: At requested position
 Models: Upper Lower Natural Shift (___ mm)
 Pictures: Packaged E-Mailed (left / right)
 Digital Scans: Please export the scans as an STL file and email it to:
 scans@respiremedical.com or submit to appropriate portal
 (Once we receive the scan you will get an e-mail confirmation)

DISCLAIMER: If the technician feels that the bite/impression and/or scans will lead to an ill fitting appliance, then he/she will contact the dental office and ask if they wish to proceed. If the Dentist wishes to still go ahead, then they are liable for the cost of a remake appliance – the laboratory is not liable.

Step 6 Check for Rush Service

CALL US to confirm your rush case PRIOR to sending in case materials.

Please Note: All devices made with traditional impressions have a 15 days turnaround time from the date we receive it! Anything before will be considered a RUSH!


- 2 Day Rush \$175
- 1 Week Rush \$125
- 1.5 Week Rush \$100
- 2 Week Rush \$75

***Please Note: EF devices can not be rushed!**

Step 7 Give Us Special Instructions

[Retention Scale]
Please check.

1 2 3
passive standard extra retentive



Signature: _____ License: _____ Date: _____